

## TECHNICAL DATA SHEET



**Enriching Biotechnology**

Telephone: +86 021 55809378

E-mail address: [marketing@bio-enriching.com](mailto:marketing@bio-enriching.com)

---

<b>Product code:</b>	P61
<b>Product name:</b>	MAG25K/Con A
<b>Surface:</b>	Con A
<b>Average Size:</b>	25µm
<b>Matrix:</b>	Cross-linked agarose
<b>Solid content:</b>	10% (v/v) medium slurry
<b>Shape:</b>	Spherical
<b>Magnetic property:</b>	Superparamagnetic
<b>Saturation magnetization:</b>	>65emu/g
<b>Magnetic separation speed:</b>	< 20s
<b>Stable in:</b>	Aqueous buffers, pH 2~12
<b>Not Stable in:</b>	Organic solvents; aqueous buffers, pH >12
<b>Product form:</b>	Suspension in PBS, pH 7.4, 0.1%BSA,0.02%NaN <sub>3</sub>
<b>Surface functional group density:</b>	~15mg Con A/mL (100%,v/v)
<b>Colour:</b>	Dark
<b>Additional remarks:</b>	Storage at 2~8 °C. Do not freeze !
<b>Expiry Date</b>	Two years after production date
<b>Protein binding capacity:</b>	10~20mg Thyroglobulin/mL (100%,v/v)
<b>Application:</b>	Glycoprotein purification