

## TECHNICAL DATA SHEET



**Enriching Biotechnology**

Telephone: +86 021 55809378

E-mail address: [marketing@bio-enriching.com](mailto:marketing@bio-enriching.com)

---

<b>Product code:</b>	ST03
<b>Product name:</b>	MPS300/Streptavidin C3
<b>Ligand</b>	Streptavidin
<b>Surface:</b>	Hydrophilic polymer coating
<b>Average Size:</b>	3 $\mu$ m
<b>Density:</b>	1.3 g/cm <sup>3</sup>
<b>Matrix:</b>	Polystyrene
<b>Solid content:</b>	10mg/mL
<b>Shape:</b>	Spherical
<b>Saturation magnetization:</b>	>8emu/g
<b>Magnetic property:</b>	Superparamagnetic
<b>Magnetic substance content:</b>	20%
<b>Stable in:</b>	Aqueous buffers, pH 4~10,
<b>Not Stable in:</b>	Aqueous buffers, pH <3 or pH >12
<b>Product form:</b>	Suspension in PBS, pH 7.4, 0.1%BSA,0.02%NaN <sub>3</sub>
<b>Sedimentation Rate:</b>	<0.08mm/min
<b>Colour:</b>	Brown
<b>Additional remarks:</b>	Storage at 2~8 °C. Do not freeze !
<b>Expiry Date</b>	Two years after production date
<b>Binding capacity for Biotin:</b>	>1100pmol/mg
<b>Technotes:</b>	Low sedimentation, binding biotin-conjugated molecules such as antibodies or DNA primers
<b>Application:</b>	Immunoassay, Preparation of single-stranded DNA, Biopanning, Sample preparation of proteins for mass Spectrometry, cell isolation